



MicroDay

Sicurezza e innovazione nel laboratorio microbiologico

Padova 29/05

Milano 31/05

AVIDITY
SCIENCE

L A B O R A T O R Y W A T E R
S Y S T E M S

Dr Mina Brett-Pitt

Confidentiality Disclaimer

- The information in this presentation is strictly confidential and proprietary to Avidity Science.
- No part of this presentation may be reproduced, stored, or transmitted in any form or by any means, electronic, mechanical, photocopying, recording, or otherwise without the prior written permission of Avidity Science.
- Other brand and product names are trademarks or registered trademarks of their respective holders.

Who we are

Avidity Science has provided water purification solutions to the research market for a combined 100+ years*



Edstrom

Started in 1969, Edstrom is the world leader in automatic watering solutions. We support animal-based research in academic, pharmaceutical and clinical segments.

triple red

Started in 1998 and acquired by Edstrom in September 2017, Triple Red has become a leading player in the U.K. by providing water purification systems, laboratory equipment, consumables and unmatched after-sales support for life sciences and healthcare customers.



Started in 1987 and acquired by Avidity in September 2018, CTC has become a leading player in the design and manufacture of filtration cartridge solutions.

We have built our company of 280+ employees through the unique combination of

3 CORE COMPETENCIES

WATER PURIFICATION AND DELIVERY

ISO 9001 2015 certified manufacturing in the USA & UK

MONITORING OF CRITICAL ENVIRONMENTS

Proprietary platforms to control watering systems and monitor research equipment

SERVICE


Nearly 25% of the company's current employees are field-based service technicians

**Edstrom and Triple Red were unified under the Avidity Science name on August 1, 2018*




Part 1: What will be covered?

- Water contaminants and quality measurements
- Laboratory Water – Classification
- Avidity Science Laboratory Water systems:
 - Puro
 - Geno
 - Pico
 - Alto
 - Duo

Global Water Cycle



Effect of contaminated water on Scientific Applications



Contaminants

1. Inorganics
2. Organics
3. Particulates\ suspended solids
4. Colloids
5. Bacteria\Viruses
6. Endotoxin
7. Dissolved gases


1. Inorganics

- Major impurity in water- simple molecules
- Mainly from water being in contact with rock salts/minerals (metals)
- Any compound that lacks a Carbon atom (CO-CO₂)
- Cations (+ve) and anions (-ve)



Inorganic water quality measurement

Conductivity $\mu\text{S-cm}$



Resistivity $\text{M}\Omega\text{-cm}$

- Meg Ohms ($\text{M}\Omega\text{-cm}$) relate to the resistivity of the water (ability to resist passage of electrical current)
- Pure water has a very high resistivity of $18.2\text{M}\Omega\text{-cm}$
- μS and $\text{M}\Omega$ have an inversely proportional relationship to each other
 $1/18.2 \text{ M}\Omega = 0.055 \mu\text{S}$ $1/100 \mu\text{S} = 0.01 \text{ M}\Omega$

Examples of U.K water quality

| | $\mu\text{S}/\text{cm}$ |
|------------|-------------------------|
| Cambridge | 715 |
| Leeds | 285 |
| Birmingham | 115 |
| Nottingham | 540 |
| Manchester | 115 |
| London | 715 |
| Oxford | 645 |
| Glasgow | 200 |

2. Organics


- Typically of biological origins- Decay of vegetal material
- Carbon containing compounds
- Total Oxidizable Organic Carbon TOOC\ TOC (PPM-PPB-PPT)
- Dissolved organics can support the growth of a wide range of microorganisms therefore disrupting several biological applications

Examples of U.K water quality

| | µS/cm | Organics |
|------------|--------------|-----------------|
| Cambridge | 715 | Low |
| Leeds | 285 | High |
| Birmingham | 115 | High |
| Nottingham | 540 | High |
| Manchester | 115 | High |
| London | 715 | Low |
| Oxford | 645 | High |
| Glasgow | 200 | Very high |



3. Particulates\ suspended solids/ 4. Colloids

- Suspended matter including hard particles (sand, rock, silt), soft particles (vegetal debris) and colloids (organic or inorganic).
- SDI: Slit Density Index
- Impact on system: Foul RO membranes interfere with valves. Analytical columns can be blocked- water turbidity



5. Bacteria\Viruses


- Chlorination does remove harmful bacteria from drinking water but it will still contain live microorganisms.



7. Dissolved gasses

- Oxygen and Nitrogen in feed water <10ppm and traces of noble gases
- CO₂ behaves differently and when it dissolves in water, it reacts to set up a series of equilibria
- $\text{CO}_2 + \text{H}_2\text{O} \leftrightarrow \text{H}_2\text{CO}_3 \leftrightarrow \text{H}^+ + \text{HCO}_3^- \leftrightarrow 2\text{H}^+ + \text{CO}_3^{2-}$
- CO₂ will pass unaltered through filters and UV treatment
- CO₂ is removed from water that has been purified by ion-exchange
- As soon as this water is dispensed into a vessel in the laboratory it comes into contact with air and CO₂ absorbed. hydrogen, bicarbonate and carbonate ions will be produced and the resistivity will fall to about 1.3 MΩ.cm

7. Dissolved gasses





L A B O R A T O R Y W A T E R –
C L A S S I F I C A T I O N

Laboratory Water - Classification

- Classified into types – I, II and III
 - Type I: Ultra-pure
 - Type II: General Lab Grade
 - Type III: Primary Grade

Laboratory Water – Type I (ASTM)

- Inorganics: 18.2 Meg Ohm-cm Resistivity
- Organics: <50ppb Total Organic Carbon
- Bacterial: <10CFU/1000ml Heterotrophic Count
- Endotoxins: <0.03 EU/ml

Laboratory Water – Type II (ASTM)*

- Inorganics: >1 Meg Ohm-cm Resistivity
- Organics: <50ppb Total Organic Carbon
- Bacterial: <10CFU/1000ml Heterotrophic Count
- Endotoxins: <0.25 EU/ml

Laboratory Water – Type III

- General product offering based on Reverse Osmosis
- Inorganics: Up to 98% rejection
- Organics: >99% rejection
- Bacterial: >99% rejection
- Particles: >99% rejection

AVIDITY
SCIENCE


W A T E R P U R I F I C A T I O N
T E C H N I Q U E S

LAB EQUIPMENT & SUPPLIES | WATER SYSTEMS | CONTROL & MONITORING | SERVICE

Water Purification Techniques


- Pre-treatment
- Reverse Osmosis (RO)
- Deionisation (DI)
- Ultra-Violet (UV)
- Ultrafiltration (UF)

Pre-treatment



Removes organics and Cl_2 by Adsorption
High Surface area


Reverse Osmosis




Very efficient and cost-effective technology for the removal of majority of impurities

The result is that **concentrate** is retained on the pressurised side of the membrane and pure **permeate** passes through to reservoir

Mixed bed deionisation




Ultra-violet (UV) irradiation



UV penetrating the cell wall of a microorganism



- Dual wavelength 185/254nm
- 254nm has bactericidal action, damaging DNA and RNA polymerase and preventing replication
- 185nm is most effective at oxidising organics, breaking large organic molecules into smaller ionised components



Ultrafiltration/Microfiltration

- An ultrafiltration (UF) filter functions as a molecular sieve
- Separates dissolved molecules according to size, size by passing the water through a very fine filter, typically with pore sizes 1-10nm
- An excellent technology for ensuring consistent ultra pure water quality with respect to particles, bacteria, pyrogens and nucleases

Avidity Science Solution

| REVERSE OSMOSIS | DEIONIZED | ULTRAPURE |
|--|--|---|
| PRIMARY GRADE WATER | GENERAL LAB GRADE WATER | 18.2 MΩ-cm |
|  <p>Puro Production Rate: 10, 20, 50 or 80L/hr. Reservoir Required</p> |  <p>Geno Production Rate: 10, 20 or 50L/hr. Reservoir Required</p> |  <p>Alto Dispensing Rate: 2L/min. Volumetric Dispense UF Options Available</p> |
| <p>Pico RO or DI Water 35L Integrated Tank Production Rate: 10 or 20L/hr.</p>  | | <p>Duo Dual Water System Production Rate: 10 or 20L/hr. Dispensing Rate: 2L/min. Volumetric Dispense UF Options Available</p>  |
|  <p>RESERVOIRS</p> |  <p>CONSUMABLES</p> |  <p>SERVICE</p> |

Puro- It speaks volumes

- Type 3 Primary Grade **RO** Water
- 10, 20, 50 or 80 litres per hour
- Choice of storage reservoir
- One easy change pre-filter
- Very low running costs
- Integral boost pump - simplifies installation and avoids unexpected cost of reduced flow rates
- Integrated leak detector
- Icon driven touchscreen – all parameters displayed



Puro



Type 3 water applications



Reservoirs

- 30, 60 and 100 litre reservoirs
- Made from Polyethylene
- 100 litre with optional pump
- Easy sanitisation – crevice free and full draining
- Wall mounting brackets for 30 and 60 litre
- Hydrophobic 0.2 μ m filter or CO₂ option



Geno. Beyond expectation

- Type 2 **General** Laboratory Grade Water
- 10, 20 or 50 litres per hour
- RO and DI technology (same RO as Puro)
- Cartridge life extended by water recirculation
- Wall, bench or reservoir mounted
- Choice of storage reservoir
- Integrated leak detector




Geno




Type 2 water applications




Normal prostate




Prostate carcinoma




Normal kidney





Renal cell carcinoma



Normal pancreas



Pancreatic carcinoma



Pico. Small in size. Big in benefits

- No other system on the market like PICO
- Type 3 or Type 2 water
- Wall, bench or reservoir mounted
- Integrated reservoir
- Very simple to install and service
- Very low running costs
- Energy consumption of <math><40\text{w}</math> when processing water
- Quite, internal boost pump and patented anti-vibration mountings



Alto - It will work around you

- Type 1 Ultrapure Water polisher
- Pre-purified water feed (tank or ring)
- On system dispense, remote dispense or both (2m tubing length)
- Wall, bench or under bench mounted
- No reservoir required
- Fast volumetric and variable dispense of Type 1 water on demand (up to 2l/min)
- Endotoxin and nuclease free (UF model)
- Dual wavelength 185/254 nm UV






Type 1 water applications




The Whole Genome Sequencing (WGS) Process

WGS is a laboratory procedure that determines the order of bases in the genome of an organism in one process. WGS provides a very precise DNA fingerprint that can help link cases to one another allowing an outbreak to be detected and solved sooner.



Duo- The best of both worlds

- Unique design:
 - Type 1 ultrapure water from dispenser
 - Type 2 water from the reservoir
- 10 or 20 l/hr Type 2 make up rates
- Space saving
- Novel tap to Type 1 water
- Fast variable dispense up to 2 l/min Type 1
- Wall, bench, under bench or reservoir mounted
- Choice of dispense – remote/system mounted
- Choice of storage reservoir



Duo












Software



Consumables

Annual low cost Endure consumable packs
Incentivise service contract with enhanced discount on consumable pack

ENVIRONMENTALLY FRIENDLY

Consumable life in Puro, Geno, Duo and Alto is optimised by utilising the latest microprocessor controls, ensuring every last drop of water is produced before they require changing. This saves the environment and makes maximum use of your budget.



Avidity Science Global Reference Sites (EMEA & Asia)



Key Life Science Accounts (LS)

| | |
|---------------------------------|--|
| UK | Phamaron UK Ltd |
| University of Cambridge | Syngenta Ltd |
| University College London (UCL) | Tecniplast UK Ltd |
| Kings College London | Immunocore Ltd |
| University of Birmingham | MRC Clinical Sciences Centre |
| Trinity College Dublin | ALS Food & Pharmaceutical |
| University of Oxford | Coca-Cola European Partners GB LTD |
| University of Southampton | Charles River Laboratories Edinburgh Ltd |
| University of St Andrews | Astex Therapeutics Ltd |
| Imperial College London | Anglian Water Services Ltd |
| AstraZeneca UK Ltd | Ireland |
| University of Manchester | France |
| University of Edinburgh | Italy |
| Cancer Research UK | Switzerland |
| Imperial College London | Germany |
| University of Hertfordshire | South Korea |
| Nuffield Healthcare Ltd | China |