Analytical Control De Mori - MicroDay 2019 - Padova & Milano



Ultra fast microbiology testing solution



What is Metis technology about?

 A rapid method for microbiological testing using flow cytometry for Quality Control in food and cosmetic industries



Key advantages




Ultimate benefit for your company



Enhance your profitability



Save time ... and money



Total gain in the production cycle: up to 12 days

Reduce leadtime

Release your product faster with maximum reliability

Decrease up to 20%
Of the storage area



A 3 steps easy to achive workflow







1

Enrichment

Minimum 12 h Optimum 24 h 2

Labeling with

METIS Fluo4Life®

viability markers

+

counterstaining with

METISol Blue®

Inhibition of natural autofluorescence

45 min

3

Counting of the viable cells with Becton Dickinson Accuri C6 PLUS ® flow-cytometer

30 s per sample



Results on **METISoft** ®: positive and negative samples at a glance.




Principle




METIS has developped proprietary fluorescent mixes for viable micro-organisms labelling:

Bacteria, Yeast, Molds









Physical parameters
Measure &
Fluorescence detection
for individual viable
microorganisms counting



Direct correlation with traditional plate count on more than hundred species






High sensitive results



PRESENCE / ABSENCE RESULTS in 24 to 48 h



Total Viable Count Sterility testing 10 to 100 times more sensitive than ATP metry



TYPE OF CONTAMINATION



a unique multiparametric detection allowing the discrimination between bacteria and yeasts & mold spores



Methods comparison





Sensitivity *

Conventionnal Plate Count

Flow Cytometry

ATP metry

pH metry

Resazurin





Flexibility



Scalable automation solutions with 2 separate devices for a minimum operator time

★ Semi-automated

Manual labeling + automatic reading Using BD Accuri C6 PLUS®





★ ★ Labelling automation high throughput

Automated labeling on an Integra pipeting robot + automatic reading using BD Accuri C6 PLUS®



★ ★★ Full automation high throughput

Automated dilution, pooling and labeling on a pipeting Hamilton workstation + automatic reading using BD Accuri C6 PLUS®





Flexibility



METIS adapts to the lab's testing throughput and can be scaled-up

564 tests in one shift 1H 2H 3H 4H 5H 6H 7H Microplate 1 labeling 94 tests +2 ctrl Reading 94 tests +2 ctrl Microplate 2 labeling 94 tests +2 ctrl Reading 94 tests +2 ctrl Microplate 3 labeling 94 tests +2 ctrl Reading 94 tests +2 ctrl Microplate 4 labeling 94 tests +2 ctrl Reading 94 tests +2 ctrl Microplate 5 labeling 94 tests +2 ctrl Reading 94 tests +2 ctrl Microplate 6 labeling 94 tests +2 ctrl Reading 94 tests +2 ctrl





A large panel of validated SOP







FOOD and BEVERAGES

Raw milk, whole milk, semi-skimmed / skimmed milk, delactosed, concentrate milk

Infant formula

Chocolate or flavoured milk (vanilla, strawberry, banana...)

Soya and vegetable drinks (rice, almond, oats...)

Dairy cream, vegetable cream

Yoghurt

Soups and high-protein content drinks (diet or medical)

Water, flavoured water and diluted fruit juices

Fruit juices with or without pulp

Carbonated beverages

Prepared baby-foods, baby-soups, compotes

Cakes, sweet-bread rolls, milk bread



COSMETIC – Personnal Care Products

Moisturing cream, anti-aging cream, solar cream, scrub, ...

Make-up: foundation, lipstick, mascara, khôl, powder, ...


Shampoo, after-shampoo, hair mask, tooth paste, dental bath,...

Oil, micellar solution, make-up removing solution, cleansing gel, ...



Reliable and user-friendly cytometer







BD Accuri

C6 PLUS®

Reliable team and know how



Over 15 years experience

R&I lab Application engineers



- Validation supports : IQ / OQ / PQ
- Detailed technical documentation : SOP/user's guide
- On-site training
- Telephone and on-site technical support
- Swapping of the flow cytometer in 24-48h



Metis key customers



Dairy







Beverage
























Cosmetic






centrepharma







Save time and money




Between 6 and 12 months Return On Investment



Take home message

Metis technology provides

- An ultra rapid microbial testing solution easy to implement
- A unique combination of speed with sensitivity
- Reach your objectives of profitability and quality together





Grazie mille



Stéphane GAVOILLE

International Sales and Marketing Director

M. +33 (0)7 89 50 72 66 sgavoille@indicia.fr

www.indicia.fr


Laetitia DA COSTA

Export Application Engineer

M. +33 (0)7 84 38 24 36 Idacosta@indicia.fr

Metis vs. Biomerieux Chemunex

- No multiparametric detection
 - → + /- results but NO discrimination between microorganisms types
- Instrument reliability :
 - → Frequent loss of laser alignment
 - → complex microfluidic



- Important volumes of sheath fluid and waste to discard
- Weaked flexibility
 - → Automated labeling and reading by the same device
 - → Testing only per batch of 25 or 50 samples in tubes : loss of reagents
 - → Max 300 samples / day
- Expensive after sale services
 - → Expensive maintenance contracts (<10 k€)
 - → "On-site intervention of service engineer" needed to repair












A propos d'Indicia

Spécialiste du contrôle de stérilité

- □ Fabricant de milieux et de produits stériles
- □ Produits et méthodes pour le contrôle de stérilité



Développement



Formulation




Fabrication




Contrôle



A propos d'Indicia













A propos d'Indicia

✓ Création

/ 3 sites

✓ Staff

✓ CA 2018

Business Model

✓ Savoir-Faire



1998

St Genis L'Argentière (Lyon) Taden (Dinan) Panazol (Limoges)



8,4 M€ (+13%)

Contract Manufacturing Organisation

Fabrication GMP de liquides stériles Microbiolologie rapide en cytométrie

Pharma/biotech, IAA, Cosmétique





