

OMNI PREP™

Programmable Multi-Sample Homogenizer

The Omni Prep™ Multi-Sample Homogenizer lets you
do the work of six people with the touch of a button!

No need to hold a motor, watch the clock, clean the probes,
or worry about sample loss or cross contamination.



- Homogenize 6 samples per cycle
- Process most samples in under one minute
- Use with Omni Tip™ reusable/disposable probes



OMNI
International
The Homogenizer Company

1000 Williams Drive, Suite 1024
Marietta, GA 30066
omni@omni-inc.com

800-776-4431
770-421-0206 Fax
www.OMNI-INC.com

OMNI PREP

Multi-Sample Homogenizer

Features

- Programmable Time and Speed (3 Stored Programs)
- Sealed or Open Vessel Processing
- Process 6 Samples Per Cycle
- Process up to 250 samples per hour
- Process Most Samples In Under One Minute
- Ultra Quiet Performance
- Clear, Removable Door
- Optional Ice Tray For Frozen Samples
- Positive Airflow That Can Be Exhausted to a HEPA-Filter or Fume Hood
- Compact Design (400mm x 200mm x 300mm)



Specifications

Processing Range: .25ml to 30ml

Speed: Variable from 500 to 30,000 rpm

Power Rating: 150 W Max

Dimensions: 400mm x 200mm x 300mm

Warranty: 1 year



Movable Rack for highly efficient processing of samples ranging from .25 to 30ml.



Tube Holders are available for 2ml, 15ml, and 50ml sample tubes.

Omni Tip plastic generator probes completely eliminate cross-contamination and cleaning between samples.



Eliminate The Tedium Of Processing One Sample At A Time

Includes: Omni Prep homogenizer and instruction manual.
Movable Rack, Tube Holders, Omni Tips and tubes sold separately.



OMNI
International
The Homogenizer Company

1000 Williams Drive, Suite 1024
Marietta, GA 30066
omni@omni-inc.com

800-776-4431
770-421-0206 Fax
www.OMNI-INC.com